



The Impulsator®

The original Impulsator® is designed for functional durability for use in a home care environment. Created in the 1980s and still used effectively today, this device delivers all of the benefits of Flow Ventilation® therapy including airway clearance and lung recruitment

Impulsator® K00012-OP-D

Designed specifically for home care and outpatient clinics, The Impulsator® allows patients to benefit from lung conditioning therapy at home. Proven safe and effective for over 30 years, the Impulsator® has provided the same Flow Ventilation®, using the Phasitron®, as the Percussionaire® hospital devices. This self-contained device is durable, easy to use, and suitable for adult and pediatric patients.

Single Patient use only.

Accessories

Phasitron Kit A50095-1D





Technical Specifications

Accessories	Part No.	
Phasitron kit	A50095-1D	
Settings	Specifications	
Aerosol flow	Liquid consumption of 1cc per minute	
Pulse/interval ratio	Average ratio 1:2.5	
Operating Conditions	Specifications	
Operating range	Temp., 0° C to 49° C (32° F to 120° F), humidity 5%-95%, Elevation 0-9842 ft	
Storage and transport range	Temp., -40° C to 60° C (-40° F to 140° F), humidity < 93% non-condensing	
Power supply	110V 50Hz	
Run time	Intermittent for 15-20 minute treatments	Caution Do not exceed 30 minutes
Display/Output	Specifications	
Pulse Amplitude	Digital display, 0 to 150 cmH2O	<p>Impulsator® display</p>
Pulse frequency	60-400 pulses per minute	
Session Usage Time	Digital display	
Filters	Description	
Red line filter	Hydrophobic	
Dimensions/Weights	Specifications	
Dimensions (WxHxD)	13" x 11.7" x 8.2"	
Weight	23 lbs	
Service life	Overhaul every 3 years from date of delivery	



PERCUSSIONAIRE®
CORPORATION

130 McGhee Road, Suite 109
Sandpoint, ID 83864
800-850-7205 | 208-263-2549
www.percussionaire.com