



**PERCUSSIONNAIRE®**  
CORPORATION

# Impulsator®



**Impulsator®**  
K00012-OP D

The **Impulsator®** is a versatile Flow Ventilation® solution built to provide lung conditioning therapy, which assists in the restoration of Functional Residual Capacity (FRC) by providing ventilation and airway clearance for patients at home.

Designed specifically for home use, the Impulsator® is a Flow Ventilation® device that is safe and provides Lung Conditioning Therapy for a wide range of home care patients.

## Accessories

**Phasitron Kit**  
A50095-1D





## Technical Specifications

Accessories	Part No.
Phasitron kit	A50095-1D
Settings	Specifications
Gas for Aerosol flow	Liquid consumption of 1cc per minute
Pulse/interval ratio	Average ratio 1:2.5
Operating Conditions	Specifications
Operating range	Temp., 5° C to 32° C (41° F to 90° F), humidity 5%-95%, Elevation 0-9842 ft
Storage and transport range	Temp., 10° C to 25° C (50° F to 77° F), humidity < 93% non-condensing
Power supply	110V 50Hz
Run time	Intermittent for 15-20 minute treatments
Display/Output	Specifications
Airway pressure	Digital display, 0 to 120 cmH2O
Pulse frequency	Digital display, 50-999 cycles/minute
Mean airway pressure	Digital display
Percussion frequency	60-400 pulses per minute
Filters	Description
Red line filter	Monitoring line filter - Hydrophobic
Dimensions/Weights	Specifications
Dimensions (WxHxD)	13" x 11.7" x 8.2"
Weight	23 lbs
Service life	Overhaul every 3 years from date of delivery



**PERCUSSIONAIRE®**  
CORPORATION

130 McGhee Road, Suite 109  
Sandpoint, ID 83864  
800-850-7205 | 208-263-2549  
www.percussionaire.com

