



PERCUSSIONAIRE®
CORPORATION

Impulsator® HC



Impulsator® HC
K00012-HC-OP D

The Impulsator® HC is a more portable Flow Ventilation® solution built to provide lung conditioning

therapy, which assists in the restoration of Functional Residual Capacity (FRC) by providing ventilation and airway clearance for patients at home.

Designed specifically for portable home use, the Impulsator® HC is a Flow Ventilation® device that is safe and provides Lung Conditioning Therapy for a wide range of home care patients.

Accessories

Phasitron Kit
A50110-1D





Technical Specifications

Accessories	Part No.
Phasitron kit	A50110-1D
Settings	Specifications
Gas for Aerosol flow	Liquid consumption of 1cc per minute
Pulse/interval ratio	Average ratio 1:2.5
Operating Conditions	Specifications
Operating range	Temp., 5° C to 32° C (41° F to 90° F), humidity 5%-95%, Elevation 0-9842 ft
Storage and transport range	Temp., 10° C to 25° C (50° F to 77° F), humidity < 93% non-condensing
Power supply	110V 50Hz
Run time	Intermittent for 15-20 minute treatments
Display/Output	Specifications
Airway pressure	No display
Pulse frequency	No display
Mean airway pressure	No display
Percussion frequency	100-500 pulses per minute, No display
Filters	Description
Filter	Inlet and Outlet filters
Dimensions/Weights	Specifications
Dimensions (WxHxD)	18" x 8" x 11"
Weight	15 lbs
Service life	Overhaul every 3 years from date of delivery.



PERCUSSIONAIRE®
CORPORATION

130 McGhee Road, Suite 109
Sandpoint, ID 83864
800-850-7205 | 208-263-2549
www.Percussionaire.com

