



**PERCUSSIONAIRE®**  
CORPORATION

# Impulsator® HC



**Impulsator® HC**  
K00012-HC-OP D

**The Impulsator® HC** lets active patients receive Flow Ventilation® therapy at home. A portable travel case allows you to discreetly get the same IPV® lung recruitment therapy you would receive in any major hospital.

Created for active users, the Impulsator® HC was designed to be portable. Proven safe and effective for over 10 years, the Impulsator® HC has provided the same Flow Ventilation® treatment, using the Phasitron®, as all Percussionaire® hospital devices. This self-contained device is easy to use, suitable for adult and pediatric patients, and may reduce hospital re-admissions.

Single Patient use only.

## Accessories

**Phasitron Kit**  
A50110-1D



# Impulsator® HC



## Technical Specifications

Accessories	Part No.
Phasitron kit	A50110-1D
Settings	Specifications
Aerosol flow	Liquid consumption of 0.5cc per minute
Pulse/interval ratio	Average ratio 1:2.5
Operating Conditions	Specifications
Operating range	Temp., 0° C to 49° C (32° F to 120° F), humidity 5%-95%, Elevation 0-9842 ft
Storage and transport range	Temp., -40° C to 60° C (-40° F to 140° F), humidity < 93% non-condensing
Power supply	110V 50Hz
Run time	Intermittent for 15-20 minute treatments  <b>Caution</b> Do not exceed 30 minutes
Display/Output	Specifications
Airway pressure	No display
Pulse frequency	No display
Mean airway pressure	No display
Percussion frequency	100-500 pulses per minute, No display
Filters	Description
Filter	Inlet and Outlet filters
Dimensions/Weights	Specifications
Dimensions (WxHxD)	18" x 8" x 11"
Weight	15 lbs
<b>Service life</b>	Overhaul every 3 years from date of delivery.



**PERCUSSIONAIRE®**  
CORPORATION

130 McGhee Road, Suite 109  
Sandpoint, ID 83864  
800-850-7205 | 208-263-2549  
[www.Percussionaire.com](http://www.Percussionaire.com)

LIT-PS-HC REV A