



The IPV®-HC allows active patients to receive TRUE-IPV® therapy at home.

The IPV®-HC
K00012-HC

The IPV®-HC Impulsator® is housed in a portable travel bag. Proven safe and effective for over 10 years, the IPV® HC has provided the same TRUE-IPV® therapy as all Percussionaire® hospital therapy devices, using the Phasitron® 5. This self-contained device is easy to use, suitable for adult and pediatric patients, and may reduce hospital re-admissions.¹

Accessories

P5-HC Phasitron® 5
Breathing Circuit Kit



¹ 2018 Nicolini et al, Safety and effectiveness of of HFCWO vs IPV® in patients with severe COPD

IPV[®] - HC

TRUE-IPV[®] THERAPY



Technical Specifications

Accessories	Part No.
Phasitron Kit	P5-HC [®]
Settings	Specifications
Aerosol Flow	Up to 10 LPM
Pulse/Interval Ratio	Automatic
Operating Conditions	Specifications
Operating Range	5° C to 40° C (41° F to 104° F) Humidity <93% non-condensing
Storage and Transport Range	-25° C to 60° C (-13° F to -140° F) Humidity <93% non-condensing
Power Input	120V 60Hz or 220v 50Hz
Run Time	Intermittent for 15-20 minute treatments
Display/Output	Specifications
Pulse Frequency	110-500 pulses per minute
Dimensions/Weights	Specifications
Dimensions (WxHxD)	45.72 cm x 45.72 cm x 27.94 cm (18" x 18" x 11")
Weight	6.8 kg (15 lbs)
Required Maintenance	Every 3 years.

**⚠ Caution Do not exceed 20 minutes
Allow a 10 minute cool down
before restart**



PERCUSSIONAIRE[®]

130 McGhee Road, Suite 109
Sandpoint, ID 83864
800-850-7205 | 208-263-2549
www.Percussionaire.com

P20059 Rev A